**METTLER TOLEDO**

## **InPro 48XX Series pH Electrode**

### **Instruction Manual**

cs	Návod k obsluze	2
da	Brugsanvisning	6
de	Betriebsanleitung	10
en	Instruction Manual	14
es	Manual de instrucciones	18
fr	Instructions d'utilisation	22
fi	Käyttöohje	26
hu	Használati utasítás	30
it	Istruzioni per l'uso	34
ja	取扱説明書	38
ko	사용 매뉴얼	42
nl	Gebruikershandleiding	46
pl	Instrukcja obsługi	50
pt	Manual de instruções	54
ru	Инструкция по эксплуатации	58
sv	Bruksanvisning	62
th	คู่มือคำแนะนำ	66
zh	说明书	70

# Řada InPro 48XX pH elektroda

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---

### Obsah

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---

### Contents

1	Introduction .....	15
2	Safety instructions .....	15
3	Product description.....	15
4	Installation and commissioning .....	16
5	Operation: Calibration of the electrode and pH transmitter.....	16
6	Maintenance .....	16
7	Removal of sources of malfunction .....	17
8	Disposal.....	17
9	Guarantee.....	17
10	Storage conditions.....	17

InPro and ISM are registered trademarks of  
Mettler-Toledo GmbH in Switzerland, the USA,  
the European Union and a further five countries.

## 1 Introduction

The METTLER TOLEDO InPro® 48XX pH-single rod measuring cell is a low maintenance, high-pressure resistant pH electrode with a pressure compensated double chamber system using a gel electrolyte and a reference system together with a silver ions barrier.

Please read through these operating instructions carefully before commissioning, in order to ensure trouble-free use.

## 2 Safety instructions

When cleaning or calibrating the electrode using solutions containing acids or alkalis both protective spectacles and protective gloves should be worn.

Do not touch the electrode at the plug-in contacts, since electrostatic discharges can damage the electronics.

We recommend that you only operate the electrode in combination with original parts from METTLER TOLEDO. Operation and maintenance should be carried out only by trained personnel and staff who have read and understood the operating instructions.

## 3 Product description

The inscription on each electrode contains the following information:

<b>METTLER TOLEDO</b>	Manufacturer of the electrode
<b>InPro® 48XX i/SG/yyy/zzz/ISM</b>	Type designation; XX = shaft material, membrane glass; i = ISM Digital/ SG = ancillary electrode/ yyy = rod length in mm/ zzz = type of temperature sensor/ ISM (analog)
<b>combination pH</b>	Type of electrode
pH 0 ... 14 (InPro 4800)	pH measurement range
pH 1 ... 14 (InPro 4801 / 4881)	
pH 1 ... 11 (InPro 4802)	
- 5 ... 130 °C	Temperature range for operation
0 ... 80 °C (InPro 4802)	
Order No. 52 00X XXX	Order number
 SEV14 ATEX 0168 X	Identification and number of certificate
 IECEx SEV 14.0025X	Identification and number of certificate
IS/I, II, III/1 / ABCDEFG/T6	FM tested
See instruction manual!	Instruction (Observe the operating instructions)

In addition each electrode is allocated a serial number on the red part of the plug-in head to enable identification.

## 4 Installation and commissioning

1. While unpacking check the electrode for mechanical damage. Please report any damage immediately to your METTLER TOLEDO supplier.
2. Remove the protective cap and rinse the electrode for a short time with deionized water. If the electrode has already been used the sensitive part can be cleaned with a soft toothbrush and a mild soap solution, and then rinsed for a short time with deionized water. After rinsing the electrode should only be dabbed. Any rubbing of the pH sensitive glass can increase the response time as a result of electrostatic charge.
3. Check the space behind the pH sensitive glass for air bubbles and remove any air bubbles that may be present by shaking the electrode slightly in a vertical direction.
4. Place the electrode into the armature as described in the instructions for the armature.
5. Connect the electrode and the pH transmitter with the corresponding cable (VP6, VP8, AK9). Observe the connections scheme that accompanies the cable or the transmitter.
6. Intelligent ISM® sensors such as InPro 480X (ISM) and InPro 48XXi enable "Plug and Measure" and expanded diagnostics. For installation, commissioning and operation of the ISM system please refer to instructions for the transmitter, module and cable.

## 5 Operation: Calibration of the electrode and pH transmitter

Before a calibration immerse the electrode for 10 minutes in a pH 7.00 or pH 4.01 buffer solution whilst connected it to the transmitter. A 2-point calibration is recommended, e.g. pH 7.00 and pH 4.01 buffer solutions. For further details please consult the operating instructions for the pH transmitter.

## 6 Maintenance

1. After each working cycle the electrode tip and diaphragm have to be carefully rinsed with deionized water. Any remaining of the solution being measured on these parts must be avoided at all costs!
2. If the electrode is not in use it has to be stored with the electrode tip and the diaphragm immersed in a tank of electrolyte (9823/3M KCl or FriscoLyf™ 9848).

3. If the electrode is stored in its armature the advices described under Item 2 apply, although the electrolyte has to be slightly modified for storage purposes (2 parts of buffer solution pH 9.2 to 10 parts of electrolyte) in order to avoid any corrosion of the armature parts.
4. If an electrode is inadvertently stored in a dry condition for several days it must be rinsed for several hours in a standard storage electrolyte.
5. The plug-in head should occasionally be checked for possible traces of moisture. If necessary clean it thoroughly with deionized water or alcohol and then dry it off carefully.

## **7 Removal of sources of malfunction**

With a soft toothbrush and water any residues can easily be removed from the electrode. Mild washing agents can also be used. Protein contaminations on the diaphragm can be removed with the cleaning solution (order no. 51 340 070).

**Warning:** Solutions containing hydrofluorides damage the pH sensitive part of the electrode.

## **8 Disposal**

It is the user's responsibility to dispose of the sensor in a professional manner. The sensor can contain electronic components that require disposal such that no persons or the environment are endangered.

## **9 Guarantee**

On manufacturing defects 12 months after delivery.

## **10 Storage conditions**

We recommend that the electrodes are stored at room temperature.

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03/2016 © Mettler-Toledo GmbH
- France**                      **Mettler-Toledo Analyse Industrielle S.A.S.**
- Germany**                      **Mettler-Toledo GmbH**
- Switzerland**                      **Mettler-Toledo (Schweiz) GmbH**
Mettler-Toledo GmbH, Process Analytics, Im Hackacker 15, CH-8902 Urdorf