# Quick Setup Guide Transmitter M800

---

## Multilingual Quick Setup Guide

---

cs	Průvodce rychlým nastavením	3
da	Lynvejledning	13
de	Quick-Setup-Leitfaden	23
en	Quick Setup Guide	33
es	Guía de configuración rápida	43
fi	Pika-asetusopas	53
fr	Guide de paramétrage rapide	63
hu	Gyorsbeállítási útmutató	73
it	Guida alla configurazione rapida	83
ja	クイックセットアップガイド	93
ko	빠른 설정 가이드	103
nl	Beknopte handleiding	113
pl	Podręcznik szybkiej konfiguracji	123
pt	Guia de Configuração Rápida	134
ru	Руководство по быстрой настройке	144
sv	Snabbinstallationsguide	155
th	คู่มือการตั้งค่าอย่างรวดเร็ว	165
zh	快速操作指南	175

---

METTLER TOLEDO

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## Převodník M800

---

### Obsah

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---

### Content

1	Safety and Intended Use	34
2	Technical Data and Further Information	37
3	Mounting	37
4	Electrical Connection	37
5	Menu Structure	38
6	Operation On-Site	39
7	Sensor Calibration	41
8	Maintenance	42
9	Disposal	42
10	EC-Declaration of Conformity	42

# 1 Safety and Intended Use

**Note:** The Quick Setup Guide is a brief operating instruction.

The M800 transmitter must only be installed, connected, commissioned, and maintained by qualified specialists e.g. electrical technicians in full compliance with the instructions in this Quick Setup Guide, the applicable norms and legal regulations.

The specialist must have read and understood this Quick Setup Guide and must follow the instructions it contains. If you are unclear on anything in this Quick Setup Guide, you must read the Operation Manual (supplied on CD-ROM). The Operation Manual provides detailed information on the instrument.

The M800 transmitter should be operated only by personnel familiar with the transmitter and who are qualified for such work.

## Intended Use

The M800 is a 4-wire transmitter for analytical measurements with 4 (0) to 20 mA output signal<sup>1)</sup>. The M800 is a multi-parameter transmitter and supports the measurements listed in the parameter fit guides. The 2-channel and 4-channel versions are compatible with ISM and flow sensors. The 1-channel version is compatible with analog and ISM sensors<sup>2)</sup>.

The M800 transmitter is designed for use in the process industries, in non-hazardous areas.

- 1) M800 Profinet doesn't have analog output
- 2) M800 Profinet supports only ISM sensors, details see Profinet parameter fit guide

## M800 parameter fit guide for the 1-channel version

Parameter	Analog	ISM
pH/ORP	•	•
pH/pNa	–	•
UniCond 2-e / UniCond 4-e	– / –	• / •
Conductivity 2-e / Conductivity 4-e	• / •	– / •
Amp. Dissolved Oxygen ppm / ppb / trace	• / • / • <sup>1)</sup>	• / • / • <sup>1)</sup>
Amp. Oxygen gas ppm / ppb / trace	• / • / • <sup>1)</sup>	• / • / • <sup>1)</sup>
Optical Dissolved Oxygen	–	• <sup>1)</sup>
Dissolved Carbon Dioxide (InPro5000i)	–	•
CO <sub>2</sub> hi (InPro5500i)	–	•
Turbidity	• (backscatter)	•

1) INGOLD sensors

## M800 parameter fit guide for the 2-channel and 4-channel versions

These versions are compatible with the following ISM and flow sensors.

Parameter	Water		Process	
	2-channel	4-channel	2-channel	4-channel
pH/ORP	•	•	•	•
pH/pNa	–	–	•	•
UniCond 2-e	•	•	•	•
Conductivity 4-e	•	•	•	•
Amp. Dissolved Oxygen ppm / ppb / trace	– / • / – <sup>2)</sup>	– / • / – <sup>2)</sup>	• / • / • <sup>1)</sup>	• / • / • <sup>1)</sup>
Amp. Oxygen gas ppm / ppb / trace	– / • / – <sup>2)</sup>	– / • / – <sup>2)</sup>	• / • / • <sup>1)</sup>	• / • / • <sup>1)</sup>
Optical Dissolved Oxygen	• <sup>2)</sup>	• <sup>2)</sup>	• <sup>1), 3)</sup>	• <sup>1), 3)</sup>
Dissolved Carbon Dioxide (InPro5000i)	–	–	•	•
CO <sub>2</sub> hi (InPro5500i)	–	–	• <sup>3)</sup>	• <sup>3)</sup>
TOC	•	•	–	–
Dissolved Ozone	•	•	–	–
Flow	•	•	–	–

1) INGOLD sensors

2) THORNTON sensors

3) 2-channel: An opt. dissolved sensor or a CO<sub>2</sub> hi sensor has to be connected to channel 2. 4-channel: Optical dissolved sensors and CO<sub>2</sub> hi sensors have to be connected to channel 2 and / or to channel 4.

## M800 Profinet parameter fit guide

Version Part no.	Process 1-ch 30 530 021	Process 2-ch 30 530 022
pH/ORP	•	•
pH/pNa	•	•
UniCond 2e/4e	•	•
Cond 4e	•	•
Amp. DO ppm/ppb/trace	•/•/•*	•/•/•*
Amp. O <sub>2</sub> gas ppm/ppb/trace	•/•/•*	•/•/•*
Optical DO	•**	•**
Dissolved Carbon Dioxide	•	•
CO <sub>2</sub> hi (thermal conductivity)	•**	•**
InPro 86x0i	•	•
Dissolved O <sub>3</sub>	–	–
Flow	–	–

\* INGOLD sensors

\*\* One optical DO or thermal conductivity CO<sub>2</sub> sensor can be used together with 2-channel transmitter

## 2 Technical Data and Further Information

The most important technical data such as supply voltage are given on the name plate either outside or inside the transmitter housing. For further technical data such as accuracy refer to the Operation Manual. This document, the Operation Manual and software are on the supplied CD-ROM. You can also download the documentation via Internet "www.mt.com/M800".

## 3 Mounting

The M800 transmitter is available as 1/2 DIN version.  
For installation drawings refer to the Operation Manual.

**DANGER! Mortal danger by electric shock or risk of electrical shock:** The maximum screw-in depth of the mounting holes in the housing is 12 mm (0.47 inch). Do not exceed maximum screw-in depth.

1. Mount the supplied cable glands at the housing.
2. Mount the transmitter. You have the following possibilities:
  - Panel mounting, wall mounting or pipe mounting

## 4 Electrical Connection


**DANGER! Mortal danger by electric shock:** Power off instrument during electrical connection.

**Note:** Do not supply power to the **Aout** terminals.


1. Switch off supply voltage.
2. Connect mains supply as follows:
  - 20 to 30 V DC: **N (-)** for Neutral and **L (+)** for Line
  - 100 to 240 V AC: **N** for Neutral and **L** for Line
3. Connect sensor, analog output signals, digital input signals and relay output signals according to the Operation Manual.

## 5 Menu Structure

\* For Profinet only

## 6 Operation On-Site

### A Start screen (example)

- 1 Changing between the channels, only 2-channel and 4-channel versions
- 2 1st line, standard configuration
- 3 2nd line, standard configuration
- 4 3rd line, depends on configuration
- 5 4th line, depends on configuration

### B Menu screen (example)

### C ISM Menu screen

Operating element	Description
	Enter Menu screen
	Enter Start screen
	Enter ISM menu
	Enter Favorite menu
	Enter Calibration menu
	Enter Configuration menu
	Return to Menu screen
	Enter next-lower menu level, here e.g. iMonitor, Messages or ISM Diagnostics
	Return to next-higher menu level
	<ul style="list-style-type: none"> <li>• Change between pages within one menu level</li> <li>• 2-channel and 4-channel: Change between the channels</li> </ul>

## 7 Sensor Calibration

 **Note:** For best process calibration results observe the following points. Take grab sample as close as possible to the measurement point of the sensor. Measure the sample at process temperature.

For a detailed description of the “Process”, “1-Point” and “2-Point” calibration method see Operation Manual of the M800 transmitter. Via iSense software you can calibrate the sensor with the “1-Point” or “2-Point” calibrate method. See Operation Manual of the iSense software.

As soon as the calibration is in progress no other calibration can be started.

### Sensor Calibration Menu

After every successful calibration different options are available. If “Adjust”, “SaveCal” or “Calibrate” is chosen, the message “Calibration saved successfully! Reinstall sensor” is displayed. Press “Done”.

Option	Analog sensors (1-channel version only)	ISM (digital) sensors
<b>Analog sensors:</b> <b>Save Cal</b>	Calibration values are stored in the transmitter and used for the measurement. Additionally, the calibration values are stored in the calibration data.	Calibration values are stored in the sensor and used for the measurement. Additionally, the calibration values are stored in the calibration history.
<b>ISM sensors:</b> <b>Adjust</b>		
<b>Calibrate</b>	The function “Calibrate” is not applicable for analog sensors.	Calibration values are stored in the calibration history for documentation, but not be used for the measurement. The calibration values from the last valid adjustment are further used for the measurement.
<b>Cancel</b>	Calibration values are discarded.	Calibration values are discarded.

## **8 Maintenance**

The transmitter requires no maintenance.

Clean the surfaces with a soft damp cloth and dry the surfaces with a cloth carefully.

## **9 Disposal**

Observe the applicable local or national regulations concerning the disposal of "Waste electrical and electronic equipment".

Disassemble the transmitter according to resources. Sort resources and supply them to recycling. Non-recyclable materials have to be disposed of in an environmental-friendly manner.

## **10 EC-Declaration of Conformity**

The EC-Declaration of conformity is part of the delivery.

ISM is a registered trademark of the Mettler-Toledo Group in Switzerland, Brazil, USA, China, European Union, South Korea, Russia and Singapore.

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© 08/2019 METTLER TOLEDO. All rights reserved.
**Mettler-Toledo GmbH**
**Mettler-Toledo Thornton, Inc.**